

RESTORATION

DENITAI
CATALOG

INNOVATIONS FOR RESTORATIVE DENTISTRY

For over 30 years, COLTENE has been a renowned provider of high grade restorative products, which have ever since developed into a focus product segment. Here the emphasis is on innovative and clinically proven products for aesthetic filling therapy, as well as a wide range of luting materials.

COLTENE pioneered the development of white tooth fillings to replace amalgam filling materials which were the standard in early days. As early as 1983, the micro-hybrid composite BRILLIANT was introduced, which is now an integral part of restorative dentistry.

Collaboration with leading universities and opinion leaders is an important success factor, which has repeatedly led to groundbreaking new developments. This is how the BRILLIANT-Esthetic Line filling concept enabled intraoral as well as extraoral restoration thanks to a special polymerisation furnace, which was far ahead of its time. To simplify shade selection, a "Duo Shade" system was developed for the SYNERGY composite. Thanks to the excellent self-blend properties, the entire shade spectrum is covered with just a few shades.

With MIRIS, the premium filling material from COLTENE, a completely new shade concept "Natural Layering" was implemented in collaboration with Dr. Didier Dietschi, University of Geneva. Here the dentine and enamel layers of the natural tooth are reproduced one to one in order to satisfy even the most demanding aesthetic requirements.

With the development of the first Direct Composite Veneering System, COMPONEER, composite enamel shells were presented for direct restoration of anterior teeth for the first time. Your patients can be provided with a natural aesthetic smile – in just one sitting! An attractive and successful treatment alternative for the patient and the dentist. In addition, a broad spectrum of outstanding adhesive systems, dependable luting cements, and also as core build-up materials, polymerisation lamps and much more besides round off the product portfolio and portray the company's innovative approach.




COMPANY FOUNDATIONS

- 1910 ROEKO GmbH & Co.KG, DE
- 1953 Vigodent, BR
- 1956 Whaledent Inc., US
- 1963 COLTENE AG, CH



MERGERS

- 1990 Whaledent Inc., US
- 1997 HYGENIC Dental, US
- 2001 DIATECH, CH
- 2002 ROEKO, DE
- 2007 CEI, US
- 2010 Vigodent, BR



TREATMENT AUXILIARIES

- 1955 Luna
- 1964 Parotisroll
- 1970 Cotton Pellets
- 1972 Surgitip
- 1975 Gelatamp
- 1995 Dental Dam Silicone non latex
- 2000 Flexi Dam non latex

INFECTION CONTROL

- 1988 BioSonic UC 1 Cleaner
- 1994 Steri-Quick
- 1997 Protecta Quick non latex
- 1999 Simplex-Plus Maxi
- 2001 BioSonic US 100 Scaler
- 2006 BioSonic UC 125
- 2012 Protecta Active Carbon

ROTARY INSTRUMENTS

- 1981 ML Diamonds
- 1992 Carbides
- 1994 Top Spin
- 2007 Speedster
- 2011 SwissFlex

**INDICATION
PRODUCT TYPE**

**FIELD OF APPLICATION
PROPERTIES**

PRODUCT NAME

ETCHING

Etching gel

35% phosphoric acid

Etchant Gel S

10

BONDING

Light-curing bond

One-component enamel adhesive

Margin Bond

10

Classical dentine enamel bond

A.R.T. Bond

10

Total-Etch dentine enamel bond

One Coat Bond

10

Self-Etching Bond

One Coat Self-Etching Bond

11

Self-Etching one component adhesive

One Coat 7.0

11

Chemically curing bond

Self-conditioning adhesive system

ParaBond

11

FILLING MATERIALS

Light-curing composites

Aesthetic filling material

MIRIS2

15

Universal nano composite

SYNERGY D6

17

Nanofilled flow material

SYNERGY D6 Flow

19

Nanohybrid composite

SYNERGY Nano Formula

21

Dual-curing composite

Flowable bulk composite

Fill-Up!

21

Accessories

Matrix bands

For amalgam fillings and glass ionomer cements

ROEKO matrix bands

22

For finishing the approximal space

ROEKO steel separating strips

22

Additional accessories

Auxiliary material

Dispenser for composite tips

22

Mixing block

22

Effect shades

Composite-based effect shades

Paint On Colors

22

Grinding and polishing strips

For smoothing the approximal space

ROEKO grinding and polishing strips

22

DIATECH

For contouring, finishing and high gloss polishing of filling materials

DIATECH SwissFlex Kit

23

DIATECH polishing kits

Silicone polishing system for pre-polishing and high gloss polishing of composites

DIATECH Composite Polishing Kit

23

Diamond polishing system for smoothing, pre-polishing and high gloss polishing of composites

DIATECH Composite Polishing Plus Kit

23

High performance diamond polishers for smoothing and high gloss polishing

DIATECH Composite 1step Polishing Kit

23

**INDICATION
PRODUCT TYPE**

**FIELD OF APPLICATION
PROPERTIES**

PRODUCT NAME

CURING LIGHTS

Light polymerisation	High performance curing light	S.P.E.C 3	25
	Standard LED	Coltolux LED	25

TEMPORARY RESTORATION

Crowns & Bridges	For the fabrication of temporary crowns and bridges	Cool Temp NATURAL	26
	A-silicone for temporary cementation	TempoSIL 2	27
Cavities	Dual-curing temporary filling material	DuoTEMP	26
	Chemically curing filling material	Coltosol F	27
Accessories	For easy removal of temporaries	Intraoral Separator	27

COMPOSITE ENAMEL SHELLS

Front tooth	Nano-hybrid composite enamel shells	COMPONEER	31
Class V	Nano-hybrid composite enamel shells	COMPONEER CLASS V	37
Accessories	For modelling composites	Modelling instrument MBC	39
	For preparing and finishing of COMPONEER	DIATECH COMPONEER Preparation and Finishing Kit	39

CORE BUILD-UP AND CEMENTS

Core build-up material	Dual curing material for core build-up and cement	ParaCore	42
Cements	Dual curing cement in the automix syringe	DuoCem	43
	Dual curing cement in the handmix	Duo Cement Plus	43
	Self-adhesive, dual curing composite-based cement	SoloCem	45

PARAPULPAL PINS

Pins	Titanium pin with round retention head	MAX Pin	46
	Titanium pin with shoulder stop	TMS LINK/ TMS LINK Plus	46
		TMS	49
	Colour-coded drills for the TMS family	KODEX drills	48

Etching
& Bonding

Restoration

LIGHT CURING ADHESIVE SYSTEMS

SELECTIVE-ETCH


Margin Bond

1st generation
Enamel is etched



A.R.T. Bond


3rd generation
Enamel is etched
Dentine is self-conditioned



TOTAL-ETCH

One Coat Bond


5th generation
Enamel is etched
Dentine is etched



SELF-ETCH

One Coat Self-Etching Bond

6th generation
Enamel is self-conditioned
Dentine is self-conditioned



One Coat 7.0


7th generation
Enamel is self-conditioned
Dentine is self-conditioned



CHEMICALLY-CURING ADHESIVE SYSTEMS


ParaBond

Enamel is etched with Etchant Gel or
Non-Rinse Conditioner
Dentine is self-conditioned



One Coat 7.0 Activator

Enamel is self-conditioned
Dentine is self-conditioned



Etchant Gel S

35 % phosphoric acid etching gel with precise placement for use on enamel and dentine.



Etchant Gel S Kit	7373
Syringes	2 × 2.5 ml
Application Needles	20 pcs

Etchant Gel S Refill	35 %	5 ml	7374
Application Needles	black	20 pcs	7876

Margin Bond

Light curing one component enamel adhesive for adhesive tooth restoration.



Margin Bond	7367
Light-cured enamel bond	5 ml

A.R.T. Bond

Classical light curing dentine/enamel bonding system for universal indications, high bond strength, withstands the test of time, scientifically recommended.



A.R.T. Bond Kit	7362
Primer A	5 ml
Primer B	5 ml
Bond	5 ml
Accessories	

Bond Refill	5 ml	7363
Primer A & B Refill	2 × 5 ml	7364

One Coat Bond

Total-Etch light curing dentin/enamel bonding system for universal indications, easy application due to non-drip gel consistency, solvent-free, prevents postoperative sensitivity.



One Coat Bond Kit	6001 1757
One Coat Bond	1.2 ml
Etchant Gel S (35 % phosphoric acid)	2.5 ml
Brushes black	50 pcs
Application Needles	20 pcs

One Coat Bond Refill	1.2 ml	7879	
One Coat Bond Refill	2 × 1.2 ml	7873	
Brushes black	50 pcs	7634	
Application Needles	black	20 pcs	7876

One Coat Self-Etching Bond

Formulated to provide an easy technique, One Coat Self-Etching Bond eliminates the need for an extra etching step. Primer 1 solubilizes the smear layer, without totally removing it as total etch bonds do. Bond 2 penetrates through the smear layer into the tubules to create the so called resin tags, which are responsible for its outstanding retention.



One Coat Self-Etching Bond Introkit		7630
Primer 1	5 ml	
Bond 2	5 ml	
Brushes yellow	50 pcs	
Brushes black	50 pcs	
Mixing well	1 pc	

One Coat Self-Etching Primer 1 Refill	5 ml	7635
One Coat Self-Etching Bond 2 Refill	5 ml	7636

Brushes yellow	50 pcs	7633
Brushes black	50 pcs	7634

One Coat Self-Etching Bond Single Dose			7631
Primer 1	50 pcs	each 0.1 ml	
Bond 2	50 pcs	each 0.1 ml	

One Coat 7.0

7th generation, single component, light curing, self-etching adhesive designed to etch, prime and bond all in one solution for "effortless bonding performance" on any surface!



One Coat 7.0 Intro Kit		7047
One Coat 7.0	5 ml	
Brushes black	50 pcs	
Mixing well	1 pc	
Cover, mixing well	1 pc	

One Coat 7.0 Activator Kit		7054
One Coat 7.0 Activator	3 ml	
Brushes black	50 pcs	

One Coat 7.0 Refill	5 ml	7049
Brushes black	50 pcs	7634

ParaBond

Chemically cured, self-conditioning adhesive system for enamel and dentine which consists of a non-rinse conditioner and a chemically curing adhesive (adhesives A and B).



ParaBond Adhesive		7486
Chemical-cured bond		
ParaBond Non-Rinse Conditioner	3 ml	
ParaBond Adhesive A	3 ml	
ParaBond Adhesive B	3 ml	

ParaBond Adhesive A & B Refill		7493
ParaBond Adhesive A	3 ml	
ParaBond Adhesive B	3 ml	

ParaBond Non-Rinse Conditioner Refill	3 ml	7494
---------------------------------------	------	------

Restoration

Filling Materials

MIRIS® 2

Highly aesthetic composite system based on the innovative Natural Layering Technique, matched to young, adult or older patients.

The Shade Guide allows the dentist intuitive shade selection and guarantees unrivalled aesthetics – also suitable for complicated cases.

8 dentine, 5 enamel and 4 effect shades

SYNERGY® D6

Universal composite using the intelligent Duo Shade system.

The optical self-blend properties of SYNERGY D6 allow a broad indication spectrum to be covered with just a few shades. Optimal processing properties, high polishability and reduced shrinkage, producing a perfect final result with a natural opalescent effect.

7 dentine and 2 enamel shades
VITA™-based



SYNERGY® D6 FLOW

Nanofilled flow material based on SYNERGY D6. Stable, flows under pressure. It is characterised by enhanced radiopacity and reduced shrinkage.

7 shades,
VITA™-based




Miris^{®2}

Artificial reality

- opalescent enamels
- brightness and vitality
- reduced shrinkage
- optimised handling properties


MIRIS2



Artificial reality

MIRIS2 takes you to another level to the ideal creation of ‘artificial reality’. The idea behind MIRIS2 is to enable tooth reconstruction with two basic compounds. The compounds have optical properties that are comparable with those of dentine and enamel so they produce natural results. For this purpose, the shades and translucency of natural teeth in adolescent, adult and elderly patients were measured and emulated. In addition to dentine and enamel shades, compounds with special effects have been developed in order to allow customised reconstruction, taking special aspects into account, and thus achieve unsurpassed aesthetics – even in highly complex cases. MIRIS2 is suitable for direct reconstruction of teeth in adolescent, adult and elderly patients.

MIRIS2 Kits



MIRIS2 Syringe Set	8431
13 Syringes (4 g / 2 ml each), 4 Syringes (2.3 g / 1.4 ml each), Accessories	
MIRIS2 Dentine Shade	8 pcs
Shade 0 (S0), Shade 1 (S1), Shade 2 (S2), Shade 3 (S3), Shade 4 (S4), Shade 5 (S5), Shade 6 (S6), Shade 7 (S7)	
MIRIS2 Enamel Shades	5 pcs
White Regular (WR), White Bleach (WB), Neutral Regular (NR), Neutral Trans (NT), Ivory Regular (IR)	
MIRIS2 Effect Shade	4 pcs
Blue (B), White (W), White Opaque (WO), Gold (G)	
Shade Guide	
Accessories	


MIRIS2 Tips Set	8465
78 Tips (0.25 g/0.125 ml each), 4 Syringes (2.3 g / 1.4 ml each), Accessories	
MIRIS2 Dentine Shade	6 pcs each
Shade 0 (S0), Shade 1 (S1), Shade 2 (S2), Shade 3 (S3), Shade 4 (S4), Shade 5 (S5), Shade 6 (S6), Shade 7 (S7)	
MIRIS2 Enamel Shades	6 pcs each
White Regular (WR), White Bleach (WB), Neutral Regular (NR), Neutral Trans (NT), Ivory Regular (IR)	
MIRIS2 Effect Shade	4 pcs
Blue (B), White (W), White Opaque (WO), Gold (G)	
Shade Guide	
Accessories	

Material selection according to the clinical situation

Teeth	Bleached			Adolescents			Adults			Elderly patients		
	Dentine	Effect	Enamel	Dentine	Effect	Enamel	Dentine	Effect	Enamel	Dentine	Effect	Enamel
Classes I-II	S2	-	WB	S2	-	WR	S3	-	WR/NR	S4	-	NR
Class III	S0-S1	-	WB	S1-S2	-	WR/WB	S2-S3	-	WR/NR/IR	S4-S5	-	NR/IR/NT
Class IV	S0-S1	B/W/WO	WB	S1-S2	B/W/WO	WR/WB	S2-S3	B/G/W/WO	WR/NR/IR	S4-S5	G/W/WO	NR/IR/NT
Class V	S0-S3	-	WB	S1-S3	-	WR/WB	S2-S5	-	WR/NR/IR	S4-S7	-	NR/IR/NT
Aesthetics	S0-S1	B/W/WO	WB	S1-S2	B/W/WO	WR/WB	S2-S3	B/G/W/WO	WR/NR/IR	S4-S5	B/G/W/WO	NR/IR/NT
Veneers	(S0-S1)	B/W/WO	WB	(S1-S2)	B/W/WO	WR/WB	S2-S3	B/G/W/WO	WR/NR/IR	S4-S5	B/G/W/WO	NR/IR/NT

The effect shades Gold, White and White Opaque can be mixed so it is possible to produce any desired opaque tint in all dentine shades (S0 to S7).

MIRIS2 Refills



Dentin Shade – Syringe		4 g
Shade 0 (S0)	1 pc	8409
Shade 1 (S1)	1 pc	8432
Shade 2 (S2)	1 pc	8433
Shade 3 (S3)	1 pc	8434
Shade 4 (S4)	1 pc	8435
Shade 5 (S5)	1 pc	8436
Shade 6 (S6)	1 pc	8437
Shade 7 (S7)	1 pc	8438

Enamel Shade – Syringe		4 g
White Regular (WR)	1 pc	8439
White Bleach (WB)	1 pc	8440
Neutral Regular (NR)	1 pc	8441
Neutral Trans (NT)	1 pc	8442
Ivory Regular (IR)	1 pc	8443

Effect Shade – Syringe		2.3 g
White (W)	1 pc	8445
White Opaque (WO)	1 pc	8446
Gold (G)	1 pc	8447
Blue (B)	1 pc	8403

Accessories		
MIRIS2 Shade Guide		8411
Application Needles	10 pcs	7874
Composite Dispenser		6671

Dentin Shade – Tips (0.25 g)	10 pcs	20 pcs
Shade 0 (S0)	8406	8408
Shade 1 (S1)	8472	8486
Shade 2 (S2)	8473	8487
Shade 3 (S3)	8474	8488
Shade 4 (S4)	8475	8489
Shade 5 (S5)	8476	8490
Shade 6 (S6)	8477	8491
Shade 7 (S7)	8478	8492

Enamel Shade – Tips (0.25 g)	10 pcs	20 pcs
White Regular (WR)	8479	8493
White Bleach (WB)	8480	8494
Neutral Regular (NR)	8481	8495
Neutral Trans (NT)	8482	8496
Ivory Regular (IR)	8483	8497

SYNERGY D6




High esthetic nano composite for everyday dentistry.


SYNERGY D6 nano-hybrid composite is the ideal modern filling material: low shrinkage, optimum handling properties, fast polishability and high resistance to operating lights. These properties make it ideal for the needs of the practising dentist.

The simple selection of shades is a special feature. The result is a selection with only 6 dentine and 2 enamel shades. This small number of shades with excellent optical blend-in properties is sufficient to cover a wide range of tooth shades.

- economical shade range and easy shade selection with a unique shade guide
- duo shade system – 2 shades in one
- excellent optical blend-in properties
- ideal handling properties
- natural opalescent effect
- reduced shrinkage




Restoration using SYNERGY D6 Dentine Shade A2/B2



Restoration using SYNERGY D6 Dentine Shade A2/B2 and Enamel Universal

* Vita™ is not a trademark of Coltène/Whaledent

SYNERGY D6 Intro Kits



SYNERGY D6 Intro Kit Syringes	8362
8 Syringes (4 g/2 ml each) and Accessories	
SYNERGY D6 Dentine Shades:	
White Bleach, A1/B1, A2/B2, A3/D3, A3.5/B3, C2/C3	
SYNERGY D6 Enamel Shades:	
Universal, White Opalescent	

Shade Guide
Accessories



SYNERGY D6 Intro Kit Tips	8363
80 Tips (0.25 g/0.125 ml each) and Accessories	
SYNERGY D6 Dentine Shades:	
White Bleach, A1/B1, A2/B2, A3/D3, A3.5/B3, C2/C3	
SYNERGY D6 Enamel Shades:	
Universal, White Opalescent	

Shade Guide
Accessories

UNIVERSAL NANO COMPOSITE



SYNERGY[®]D6
beauty made simple

SYNERGY D6 Refills



SYNERGY D6 Syringes Refills	4 g/2 ml	
Dentin White Bleach (WB)	1 pc	8365
Dentin A1/B1	1 pc	8366
Dentin A2/B2	1 pc	8367
Dentin A3/D3	1 pc	8368
Dentin A3.5/B3	1 pc	8369
Dentin C2/C3	1 pc	8370
Dentin A4/C4	1 pc	8384
Enamel Universal	1 pc	8371
Enamel White Opalescent (WO)	1 pc	8372

SYNERGY D6 Tips Refills	0.25 g/0.125 ml	
Dentin White Bleach (WB)	10 pcs	8375
Dentin A1/B1	10 pcs	8376
Dentin A2/B2	10 pcs	8377
Dentin A3/D3	10 pcs	8378
Dentin A3.5/B3	10 pcs	8379
Dentin C2/C3	10 pcs	8380
Dentin A4/C4	10 pcs	8383
Enamel Universal	10 pcs	8381
Enamel White Opalescent (WO)	10 pcs	8382

SYNERGY D6 Accessories		
Shade Guide	1 pc	8410
Composite Dispenser	1 pc	6671

SYNERGY D6 Flow



SYNERGY D6 Flow is a nano-filled, radiopaque flowable composite based on the reliable and successful technology of SYNERGY D6. SYNERGY D6 Flow provides improved radiopacity and reduced shrinkage. The material is compact, but flows under pressure, provides easy handling and precise application.

SYNERGY D6 Flow is the perfect complement to SYNERGY D6 and allows an ideal material selection.



SYNERGY D6 Flow Intro Kit	7650	
A2/B2, A3/D3, White Bleach (WB), Enamel Universal	4 × 2.3 g	
Bent Application Needles	32 pcs	

SYNERGY D6 Flow Dentin Intro Kit	7660	
SYNERGY D6 Flow A1/B1, A2/B2, A3/D3, A3.5/B3	4 × 2.3 g	
Bent Application Needles	32 pcs	




SYNERGY D6 Flow Refills	2.3 g	
A1/B1	1 pc	7651
A2/B2	1 pc	7652
A3/D3	1 pc	7653
A3.5/B3	1 pc	7654
A4/C4	1 pc	7655
White Bleach (WB)	1 pc	7656
Enamel Universal	1 pc	7657

SYNERGY D6 Accessories		
Shade Guide SYNERGY D6	1 pc	8410
SYNERGY D6 Flow Application Needles	20 pcs	8469

Hole in One

Filling in a single step with **Fill-Up!™**



Fill-Up!™ - in a single step to give a perfect result

- optimal depth polymerisation with minimal shrinkage due to dual curing system restoration
- guaranteed single-layer technique - even in very deep cavities of 10 mm
- optimised sealing of margins - reduced postoperative sensitivity
- universal shade in a convenient Automix syringe for efficient placement

Deep. Fast. Perfect.

info.uk@coltene.com | www.coltene.com



COLTENE

SYNERGY Nano Formula

Radiopaque, highly-filled, nano-hybrid composite available in two consistencies.

SYNERGY Nano Formula Duo Shade

Smooth, non-sticky universal composite for all applications, covering a wide spectrum of tooth shades with fewer syringes.



Duo Shade Syringes Refills		4 g
A1/D2	1 pc	7821
A2/B2	1 pc	7822
A3/D3	1 pc	7823
A3.5/B3	1 pc	7824
A4/M5	1 pc	7825
C2/C3	1 pc	7826

Duo Shade Tips Refills		0.25 g
A1/D2	20 pcs	7831
A2/B2	20 pcs	7832
A3/D3	20 pcs	7833
A3.5/B3	20 pcs	7834
A4/M5	20 pcs	7835
C2/C3	20 pcs	7836

SYNERGY Nano Formula Compact

Packable composite, an ideal aesthetic alternative to amalgam for larger restorations.



Compact Syringes Refills		4 g
A2/B2	1 pc	7491
A3/D3	1 pc	7492

Compact Tips Refills		0.25 g
A2/B2	20 pcs	7837
A3/D3	20 pcs	7838
A3.5/B3	20 pcs	7839

Fill-Up!

Radiopaque, flowable composite material with antibacterial zinc oxide used for bulk filling therapy in the posterior dental region.

Owing to its dual curing properties, fillings can be administered and cured in a single step in arbitrary thickness.



Fill-Up! Intro Kit		6001 9342
Fill-Up! Universal	2 × 4.5 g	
Etchant Gel S	1 × 2.5 ml	
ParaBond Adhesive A & B	2 × 3 ml	
Mixing Tip Short Super Fine	6 pcs	
Mixing Tip Short Fine	8 pcs	

Mixing Tip Short Super Fine		6747
Mixing Tip, Ø 1 mm	40 pcs	

Fill-Up! Economy Kit		6001 9355
Fill-Up! Universal	5 × 4.5 g	
Mixing Tip Short Super Fine	15 pcs	
Mixing Tip Short Fine	20 pcs	

Mixing Tip Short Fine		6759
Mixing Tip, Ø 1.8 mm	40 pcs	

Fill-Up! Refill		6001 9341
Fill-Up! Universal	1 × 4.5 g	
Mixing Tip Short Super Fine	3 pcs	
Mixing Tip Short Fine	4 pcs	

Etchant Gel S Kit		7373
Syringes	2 × 2.5 ml	
Application Needles	20 pcs	

ParaBond Adhesive A & B Refill		7493
ParaBond Adhesive A	3 ml	
ParaBond Adhesive B	3 ml	

ROEKO Matrix Bands

For use with amalgam fillings and glass ionomer cement.



ROEKO Steel Matrix Band soft			
Stainless steel, extra fine (0.05 mm), can be applied interproximal			
Width 5 mm	3 m	591 505	
Width 6 mm	3 m	591 506	
Width 7 mm	3 m	591 507	

ROEKO Strip		600 009
Transparent Matrix bands made from polyester for composite fillings		
Width 9 mm		15 m

ROEKO Steel Matrix Band (regular)			
Stainless steel, regular (0.03mm)			
Width 5 mm	3 m	590 305	
Width 6 mm	3 m	590 306	
Width 7 mm	3 m	590 307	

Stainless steel, regular (0.05mm)			
Width 5 mm	3 m	590 505	
Width 6 mm	3 m	590 506	
Width 7 mm	3 m	590 507	

ROEKO Steel separating strips

For approximal separation.



4 mm	12 pcs	570 004
6 mm	12 pcs	570 006
8 mm	12 pcs	570 008

Additional accessories



Dispenser for Composite Tips		6671
Mixing Pad	1 pc	7080

Paint On Colors

Composite effect shades for perfect aesthetic enhancement of every restoration. Especially well suited for characterisation of tooth necks and fissures, for replicating chalk spots and enamel cracks, as well as for colouring dentine.

Paint On Colors Kit			7330
White, Red, Grey, Yellow, Brown, Blue	6 × 3 ml		
White Opaque	2 × 3 ml		

Paint On Colors		7338
White Opaque		1 × 3 ml

ROEKO Abrasive- and Polishing Strips

For smoothing of interproximal spaces not accessible or difficult to access with other instruments. Made of abrasion resistant polyester for final finishing of the filling.



30 micron, Width 5 mm	10 m	560 030
40 micron, Width 5 mm	10 m	560 040
60 micron, Width 5 mm	10 m	560 060

DIATECH SwissFlex Kit



Colour	Indication	Description
Black	contouring	coarse (70µm)
Blue	finishing	medium (50µm)
Red	polishing	fine (30µm)
White	high lustre	ultra-fine (5µm)

For contouring, finishing and high gloss polishing of composite, amalgam, glass ionomer cement, semi-precious metal and precious metal.

DIATECH SwissFlex Kit			250095AA		
Content			Content		
3 × 702RA		232018AA			
25 × 703UM		232019AA	25 × 707UM		232023AA
25 × 704UM		232020AA	25 × 708UM		232024AA
25 × 705UM		232021AA	25 × 709UM		232025AA
25 × 706UM		232022AA	25 × 710UM		232026AA
20 × 700UM					232016AA
20 × 701UM					232017AA

DIATECH Polishing Kits



Two-step silicon polishing system for pre-polishing and high lustre polishing of composites e.g. SYNERGY D6, SYNERGY D6 FLOW or Cool Temp NATURAL.

DIATECH Composite Polishing Kit			250014AA		
Content			Content		
2 × 2101RA		230055AA	2 × 2201RA		230072AA
2 × 2103RA		230059AA	2 × 2203RA		230076AA
2 × 2111RA		230063AA	2 × 2211RA		230080AA
2 × 2127RA		230069AA	2 × 2227RA		230086AA



Two-step diamond polishing system for smoothing, pre- and high lustre polishing of composites, e.g. SYNERGY D6 or MIRIS2.

DIATECH Composite Polishing Plus Kit			250092AA		
Content			Content		
1 × 2301RA		230291AA	1 × 2401RA		230299AA
1 × 2303RA		230293AA	1 × 2403RA		230301AA
1 × 2311RA		230295AA	1 × 2411RA		230303AA
1 × 2327RA		230297AA	1 × 2427RA		230305AA

High performance polishers for smoothing and high gloss polishing of all composite materials.



DIATECH Composite 1step Polishing Kit			250093AA		
Content			Content		
2 × 9101RA		230264AA	2 × 9111RA		230267AA
2 × 9103RA		230265AA	2 × 9121RA		230271AA

LED CURING LIGHT



S.P.E.C. 3[®]
LED

Performance Without Compromise

S.P.E.C. 3 meets the highest demands on performance in the dental practice. Due to its high light intensity, the polymerisation times can be reduced efficiently depending on the indication and the material.

- multiple modes for your dental curing needs
- high intensity output for rapid polymerization
- no wasted chair time
- patient convenience
- intuitive and simple controls for ease of use

S.P.E.C 3



Performance Without Compromise

S.P.E.C 3 meets the highest demands of performance in the dental practice. The stable, lightweight aluminium handpiece and the 360° rotating light guide provide convenient handling. Three different curing modes provide performance without compromise, enabling variable use to meet the individual requirements of dentists. A high performance, interchangeable lithium polymer battery provides over 300 10s cures between recharges.

S.P.E.C 3 LED Kit Euro 6001 3942

S.P.E.C 3 LED Euro Plug	
Light Shield	20 pcs
Barrier Sleeves	100 pcs
S.P.E.C 3 Charger Base, Power Supply	

S.P.E.C 3 LED Kit UK 6001 3943

S.P.E.C 3 LED UK Plug	
Light Shield	20 pcs
Barrier Sleeves	100 pcs
S.P.E.C 3 Charger Base, Power Supply	

Light Shield 8 × 11 mm	20 pcs	6001 3948
Light Shield 11 × 11 mm	20 pcs	6001 4360
Barrier Sleeves	100 pcs	6001 3949
Light Guide, Turbo-Tip 8 mm		6001 3950
Light Guide, 11 mm		6001 3951
Battery Module		6001 3952
S.P.E.C 3 Charger Base		6001 3953
Power Supply		6001 3955
Eye Shield		6001 4444

Coltolux LED



The Coltolux LED combines a slim curing probe with high output in a cordless lightweight penstyle design. The Coltolux LED's heat sink technology reduces the amount of heat emission and eliminates the need for active cooling.

Coltolux LED Kit C7970

Handpiece	
Charger Base	
Power Supply	
Light Shields	3 pcs
Light Lenses	5 pcs
Barrier Sleeves	100 pcs

Light Lenses	25 pcs	C7974
Barrier Sleeves	1000 pcs	C7957
Light Shields	100 pcs	C7944
Light Lenses, attachable	25 pcs	C7972
Light Shields, attachable	25 pcs	C7990

Cool Temp NATURAL

Temporary crowns and bridge material based on bis-acryl composite, available in automix cartridges and automix syringes.

- 3 fluorescent shades
- low shrinkage
- excellent polishability
- low polymerisation temperature



Cool Temp Starter Kit	5800
Cool Temp cartridge A2	85 g
Cool Temp cartridge A3.5	85 g
Mixing Tip Blue pointed	20 pcs
SYNERGY D6 Flow A2/B2	2.3 g
SYNERGY D6 Flow A3.5/B3	2.3 g
Bent Flow Application Needles	4 pcs

Cool Temp Refill Cartridge A1	5815
Cool Temp Cartridge A1	85 g
Mixing Tip Blue pointed	10 pcs

Cool Temp Refill Cartridge A2	5805
Cool Temp Cartridge A2	85 g
Mixing Tip Blue pointed	10 pcs

Cool Temp Refill Cartridge A3.5	5810
Cool Temp Cartridge A3.5	85 g
Mixing Tip Blue pointed	10 pcs

Cool Temp Automix Syringe 5 ml Refill A1	5820
Cool Temp Automix Syringe 5 ml A1	2 pcs
Mixing Tip Brown pointed	8 pcs

Cool Temp Automix Syringe 5 ml Refill A2	5824
Cool Temp Automix Syringes 5 ml A2	2 pcs
Mixing Tip Brown pointed	8 pcs


Cool Temp Automix Syringe 5 ml Refill A3.5	5825
Cool Temp Automix Syringe 5 ml A3.5	2 pcs
Mixing Tip Brown pointed	8 pcs

Accessories		
Dispenser 50/75 ml	1 pc	6900
Mixing Tip Blue pointed	40 pcs	7774
Mixing Tip Brown pointed	40 pcs	6748
Separator Intraoral	5 ml	7235

DuoTEMP

Dual curing temporary filling material, eugenol-free and radiopaque for temporary fillings, temporary fillings for endodontic applications and temporary fillings for inlay- and onlay-preparations. Available in universal white.

- dual curing
- excellent marginal fit
- easy to remove
- one material for several indications



DuoTEMP Intro Kit	5830
Syringes	2 x 5 g
Separator Intraoral	5 ml

DuoTEMP Single Pack	5831
Syringe	5 g

DuoTEMP Eco Pack	7090
Syringes	5 x 5 g

Coltosol F




Temporary, eugenol-free filling material for temporary fillings and temporary root canal seals. Available in universal white.

- fluoride-releasing X-ray opaque
- easy adaption and adhesion in the cavity
- no mixing, easy to remove
- self-curing under the action of saliva



Coltosol F Single Pack	5906
Coltosol F in Jar	38 g



Coltosol F in Syringes	5930
Syringes	5 × 8 g

Coltosol F Eco Pack	5911
Coltosol F in Jar	3 × 38 g

Intraoral Separator



For intraoral use for easy removal of temporaries.




Intraoral Separator	7235
Bottle	5 ml

TempoSIL 2




Temporary PVS cement material for different surfaces and indications.

- secure adhesion of provisional restorations to natural tooth structure, core build up material and metal alloys as well as zirconia
- simple application as well as easy removal
- fast setting time
- two shades for esthetic applications



TempoSIL 2 Trial Pack	3686
TempoSIL 2 dentin	5 ml
TempoSIL 2 white	5 ml
Mixing Tip Brown long fine	15 pcs



TempoSIL 2 Dentin Refill	6728
TempoSIL 2 dentin	4 × 5 ml

TempoSIL 2 Dentin Intro Kit	6727
TempoSIL 2 dentin	2 × 5 ml
Mixing Tip Brown long fine	15 pcs

TempoSIL 2 White Refill	6730
TempoSIL 2 white	4 × 5 ml

TempoSIL 2 White Intro Kit	6729
TempoSIL 2 white	2 × 5 ml
Mixing Tip Brown long fine	15 pcs

Mixing Tip Brown long fine	6749
Mixing Tips	40 pcs

Composite Enamel Shells

Restoration



COMPONEER®

The direct composite veneering system for highly aesthetic anterior restorations.

COMPONEER simplifies build-up and surface design in the anterior region. At the same time, it allows minimally invasive preparation of the tooth substance and individual matching. The innovative micro-retentive surface increases wettability and strengthens sustainable bonding.

3 enamel shades



COMPONEER® CLASS V

The first composite veneering technique for the successful restoration of tooth necks.

As a result of the good wetting and marginal adaptation plaque accumulation, discolouration and secondary caries are minimised. This also optimises shrinking at the margin.

3 enamel shades

COMPONEER®

Surprise your patients with a new smile –
in only one session!

Innovative. Time-saving. Surprisingly easy. COMPONEER is the Direct Composite Veneering System used for quick, easy and save restorations of single or multiple teeth. This offers new perspectives for you and your patients. So both of you have a reason to smile.

www.componeer.info



THE SMILE TO GO.



Die COMPONEER® benefits:

No laboratory required | One session | Naturally aesthetic corrections using freehand technique | Easy application with prefabricated composite veneers | Brilliant result | Attractive added value

www.coltene.com

 COLTENE

COMPONEER



COMPONEER - The Smile to go.

COMPONEER Direct Composite Veneers are a whole new category of veneers. The polymerized prefabricated nano-hybrid composite enamel shells eliminate laborious configuration of anatomical shape and surface, as well as time-consuming preparation by freehand technique. At the same time, COMPONEER can be adapted as with the individual situation anytime by using composite. The extremely thin veneer thicknesses of 0.3 mm enable preparation with minimal damage to tooth substance.

The polished and naturally designed surface gives the restoration its vital character. The micro-retentive inner surface increases wettability and strengthens sustainable bonding.

Whether incisal edge extension, covering tooth discolouration or malpositioning of teeth: with COMPONEER, aesthetic corrections are easy to achieve and finish in a single session.

Properties in clinical use.

1. Initial situation: discoloured teeth, extensive caries, malposition, incorrect axial alignment.



2. Select enamel and dentine shade with shade guide prior to isolation, ideally under daylight lamp.



3. The perfectly matching COMPONEER is selected using the bluish transparent Contour Guides.



4. Rubberdam adhesion technique facilitates functional working.



5. Prior to opening the sealed COMPONEER, the shape is verified again following preparation.



6. Careful grinding of the COMPONEER to give an exact fit.



7. The composite is distributed on the tooth using the MBC modelling instrument. The placer is used for stress-free application of the COMPONEER to its final position.



8. Final check for perfect edges before removing the rubberdam.



COMPONEER SYSTEM KITS



1301

COMPONEER PREMIUM System Kit (84 COMPONEER)	Tips 1300	Syringes 1301
each upper Set: 11, 12, 13, 21, 22, 23		
1 × UJ S	Enamel Universal	
1 × UJ M	Enamel Universal	
2 × UJ L	Enamel Universal	
1 × UJ XL	Enamel Universal	
1 × UJ S	Enamel White Opalescent	
1 × UJ M	Enamel White Opalescent	
2 × UJ L	Enamel White Opalescent	
1 × UJ XL	Enamel White Opalescent	
each lower Set: 31, 32, 33, 41, 42, 43		
1 × LJ S	Enamel Universal	
1 × LJ M	Enamel Universal	
1 × LJ S	Enamel White Opalescent	
1 × LJ M	Enamel White Opalescent	
Material and Accessories*		

* Material and Accessories
COMPONEER Contour Guide, each size (upper: S, M, L, XL; lower: S, M)
One Coat Bond (1 × 1.2 ml), Etchant Gel S (1 × 2.5 ml)
SYNERGY D6 Dentin A1/B1, A2/B2; SYNERGY D6 Enamel White Opalescent, Enamel Universal SYNERGY D6 Flow Enamel Universal (1 × 2.3 g Syringe) (each colour: 10 × 0,25 g Tips; 1 × 4 g Syringe)
1 Shade Guide Componeer-SYNERGY D6,
15 Bent Appl. Needles "Flow"
1 MBC Modelling Instrument
1 Holder und 10 Holder Caps
1 Placer, 2 Placer Adapter, 10 Placer Stamps
1 Brushholder black and 50 Brushes
1 Comprepol Plus Cup and Point, Composhine Plus Cup and Point (1 pc each)
1 SwissFlex Trial Kit (1 Mandrel, 11 Discs)



1305

COMPONEER BASIC System Kit (36 COMPONEER)	Tips 1304	Syringes 1305
each upper Set: 11, 12, 13, 21, 22, 23		
1 × UJ M	Enamel Universal	
1 × UJ L	Enamel Universal	
1 × UJ M	Enamel White Opalescent	
1 × UJ L	Enamel White Opalescent	
each lower Set: 31, 32, 33, 41, 42, 43		
1 × LJ M	Enamel Universal	
1 × LJ M	Enamel White Opalescent	
Material and Accessories*		



6001 4222


COMPONEER INTRO System Kit (24 COMPONEER)	Tips 6001 4221	Syringes 6001 4222
each upper Set: 11, 12, 13, 21, 22, 23		
1 × UJ M	Enamel Universal	
1 × UJ L	Enamel Universal	
1 × UJ M	Enamel White Opalescent	
1 × UJ L	Enamel White Opalescent	
Material and Accessories*		

COMPONEER Refills UJ

A wide selection of sizes and shades

Choose from the wide range of maxillary veneers (maxilla) available in sizes S, M, L and XL, as well as mandibular veneers (mandible) in sizes S or M. Three natural shades allow maximum aesthetics for every indication.

Upper jaw						
FDI _i	13	12	11	21	22	23



Upper jaw		
FDI _i	14	24



COMPONEER Sets UJ	Dentin Bleach Opaque	Enamel White Opalescent	Enamel Universal
COMPONEER Set UJ S , 6 Veneers: 11, 12, 13, 21, 22, 23		1311	1316
COMPONEER Set UJ M , 6 Veneers: 11, 12, 13, 21, 22, 23	6001 4500	1312	1317
COMPONEER Set UJ L , 6 Veneers: 11, 12, 13, 21, 22, 23	6001 4509	1313	1318
COMPONEER Set UJ XL , 6 Veneers: 11, 12, 13, 21, 22, 23	6001 4516	6001 4010	6001 4012
COMPONEER Set UJ Premolar S , 2 Veneers: 14/24	6001 4523	6001 4393	6001 4392
COMPONEER Set UJ Premolar L , 2 Veneers: 14/24	6001 4524	6001 4396	6001 4394

COMPONEER Refills UJ			
COMPONEER Refill UJ S , 2 Veneers	11		1321 1351
	12		1322 1352
	13		1323 1353
	21		1324 1354
	22		1325 1355
	23		1326 1356
COMPONEER Refill UJ M , 2 Veneers	11	6001 4502	1327 1357
	12	6001 4504	1328 1358
	13	6001 4505	1329 1359
	21	6001 4506	1330 1360
	22	6001 4507	1331 1361
	23	6001 4508	1332 1362
COMPONEER Refill UJ L , 2 Veneers	11	6001 4510	1333 1363
	12	6001 4511	1334 1364
	13	6001 4512	1335 1365
	21	6001 4513	1336 1366
	22	6001 4514	1337 1367
	23	6001 4515	1338 1368
COMPONEER Refill UJ XL , 2 Veneers	11	6001 4517	6001 3998 6001 4004
	12	6001 4518	6001 3999 6001 4005
	13	6001 4519	6001 4000 6001 4006
	21	6001 4520	6001 4001 6001 4007
	22	6001 4521	6001 4002 6001 4008
	23	60014522	6001 4003 6001 4009

COMPONEER Refills LJ

Lower jaw						
FDI ₁	43	42	41	31	32	33



COMPONEER Sets LJ	Dentin Bleach Opaque	Enamel White Opalescent	Enamel Universal
COMPONEER Set LJ S , 6 Veneers: 31, 32, 33, 41, 42, 43		1314	1319
COMPONEER Set LJ M , 6 Veneers: 31, 32, 33, 41, 42, 43		1315	1320

COMPONEER Refills LJ			
COMPONEER Refill LJ S , 2 Veneers	31	1339	1369
	32	1340	1370
	33	1341	1371
	41	1342	1372
	42	1343	1373
	43	1344	1374
COMPONEER Refill LJ M , 2 Veneers	31	1345	1375
	32	1346	1376
	33	1347	1377
	41	1348	1378
	42	1349	1379
	43	1350	1380



Order a COMPONEER Kit and receive your patient info kit free of charge. The order card is in your COMPONEER kit.

COMPONEER Info Kit D	6001 3798
DVD and CD, Patient brochure, Poster	



COMPONEER Accessories



Additional material for a perfect outcome

Material	
Etchant Gel S	see page 10
One Coat Bond	see page 10
SYNERGY D6	see page 17-19
SYNERGY D6 Flow	see page 19

Accessories		
Instruments		
MBC Modelling Instrument	1 pc	6001 1110
Holder	1 pc	6001 1105
Holder Caps	40 pcs	6001 1106
Placer	1 pc	6001 1107
Placer Adapter	5 pcs	6001 1108
Placer Stamps	50 pcs	6001 1109
Brushholder black	1 pc	6001 1111
Brushes	500 pcs	6001 1112

Shade Guides		
Shade Guide Synergy D6	1 pc	8410
Shade Guide COMPONEER-Synergy D6 BL/O	1 pc	6001 3812

Contour Guides	
UJ Small	7000
UJ Medium	7001
UJ Large	7002
UJ X-Large	6001 4011
UJ Premolar	6001 4401
LJ Small	7003
LJ Medium	7004

Polishing		
DIATECH SwissFlex Kit	see pag. 23	250 095AA
DIATECH Composite Polishing Kit	see pag. 23	250 014AA
DIATECH Composite Polishing Plus Kit	see pag. 23	250 092AA
DIATECH Composite 1step Polishing Kit	see pag. 23	250 093AA

COMPONEER Preparation and Finishing Kit		6001 3844
Preparation, Finishing, Polishing		
1 × G898-314-016-11-ML	200736AA	1 × G860-314-012-10-UF 60013866
1 × G368-314-023-5-ML	200015AA	1 × 2311RA 230295AA
1 × G368-314-023-5-F	200014AA	1 × 2303RA 230293AA
1 × G864-314-014-12-F	200572AA	1 × 2411RA 230303AA
1 × G801-314-023-F	200095AA	1 × 2403RA 230301AA
1 × G863-314-012-10-UF	200561AA	1 × 7003RA 232006AA

Veni.



Vidi.

COMPONEER® CLASS V

The first Composite Veneering System for the successful restoration of tooth necks

- easier and safer application
- optimal wetting of the surface
- improved marginal adaptation
- long-lasting restoration



COMPONEER CLASS V





COMPONEER CLASS V – The unique way of restoring tooth necks.

Experience has shown that defects of the tooth neck can often only be restored insufficiently using the classical methods.

COMPONEER CLASS V, the first composite veneering technique for the successful restoration of tooth necks, represents a new and innovative solution.

The specially treated nano-hybrid composite tooth neck veneers allow easy and safe restoration of defects. The special positioning guide on the veneer enables controlled and reliable application. The positioning guide can be conveniently removed later at the predetermined breaking point. The evenly applied pressure created when placing the tooth neck shell improves wetting between tooth and composite. This renders good marginal adaptation to minimise accumulation of plaque, discolouration and secondary caries. This also optimises shrinking at the margin.

Properties in clinical handling



COMPONEER CLASS V System Kit



6001 8046



6001 8056

COMPONEER CLASS V Intro Kit	6001 8040
20 COMPONEER CLASS V	
Shades: Dentine A2/B2, A3/D3	5 pcs each S and L
One Coat Self-Etching Bond, Primer 1	5 ml
One Coat Self-Etching Bond, Bond 2	5 ml
Brushes yellow	50 pcs
Brushes black	50 pcs
Mixing well	1 pc
COMPONEER CLASS V Shade Guide	
COMPONEER CLASS V Contour Guides S, S/L, L, L/L	

COMPONEER CLASS V Premium Kit	6001 8056
40 COMPONEER CLASS V	
Shades: Dentine A1/B1, A2/B2, A3/D3, A4/C4, Enamel Universal	2 pcs each S, S/L, L, L/L
One Coat Self-Etching Bond, Primer 1	5 ml
One Coat Self-Etching Bond, Bond 2	5 ml
Brushes yellow	50 pcs
Brushes black	50 pcs
Mixing well	1 pc
COMPONEER CLASS V Shade Guide	
COMPONEER CLASS V Contour Guides S, S/L, L, L/L	

COMPONEER CLASS V Refills



COMPONEER CLASS V Refills		
4 Veneers A1/B1	S	6001 8009
4 Veneers A2/B2	S	6001 8010
4 Veneers A3/D3	S	6001 8011
4 Veneers A4/C4	S	6001 8012
4 Veneers UNIV.	S	6001 8013
4 Veneers A1/B1	S/L	6001 8015
4 Veneers A2/B2	S/L	6001 8016
4 Veneers A3/D3	S/L	6001 8017
4 Veneers A4/C4	S/L	6001 8018
4 Veneers UNIV.	S/L	6001 8019
4 Veneers A1/B1	L	6001 8020
4 Veneers A2/B2	L	6001 8021
4 Veneers A3/D3	L	6001 8022
4 Veneers A4/C4	L	6001 8023
4 Veneers UNIV.	L	6001 8024
4 Veneers A1/B1	L/L	6001 8025
4 Veneers A2/B2	L/L	6001 8026
4 Veneers A3/D3	L/L	6001 8027
4 Veneers A4/C4	L/L	6001 8028
4 Veneers UNIV.	L/L	6001 8029

COMPONEER CLASS V Accessories

Additional material for a perfect outcome



6001 1110



6001 3844

Material		
One Coat Self-Etching Bond, Primer 1, Refill	5 ml	7635
One Coat Self-Etching Bond, Bond 2, Refill	5 ml	7636
Brushes yellow	50 pcs	7633
Brushes black	50 pcs	7634
SYNERGY D6	see pages 17-19	
SYNERGY D6 Flow	see page 19	

Accessories		
Instruments		
MBC Modelling Instrument	1 pc	6001 1110

Shade Guides / Contour Guides		
COMPONEER CLASS V Shade Guide		6001 8250
COMPONEER CLASS V Contour Guides Set	1 pc each S, S/L, L, L/L	6001 8422

Preparation, Finishing, Polishing		
COMPONEER Preparation and Finishing Kit	see page 35	6001 3844

Polishing		
DIATECH SwissFlex Kit	see page 23	250095AA
DIATECH Composite Polishing Kit	see page 23	250014AA
DIATECH Composite Polishing Plus Kit	see page 23	250092AA
DIATECH Composite 1step Polishing Kit	see page 23	250093AA

Core Build-up
Cements
Restoration
Parapulpal pins



	SoloCem®	DuoCem™	ParaCore®
Description	Self-adhesive, dual-curing composite-based cement	Dual-curing universal composite-based cement	Dual-curing cement for fixation and core build-ups based on composite
Adhesive system	No adhesive necessary	ParaBond A.R.T. Bond One Coat Bond	ParaBond A.R.T. Bond One Coat Bond
Shades	Dentine White opaque	Dentine White opaque	Dentine White opaque Transparent

INDICATIONS

Crowns & bridges	+++	+++	++
Veneer		+*	
Adhesive bridges		+++*	
Inlay / Onlay	+++	+++	++
Posts	+++	+++	+++
Implant abutments	++	+	
Core build-up			+++


- + possible
- ++ recommended
- +++ highly recommended

* Translucent materials require light-curing adhesive systems like A.R.T. Bond or One Coat Bond.

ParaCore



Dual curing core material and resin cement in 5 ml automix-syringe.



ParaCore 5ml Intro Kit	5885
ParaCore Dentin	1 × 5 ml
ParaCore White	1 × 5 ml
ParaCore Dentin SLOW	1 × 5 ml
ParaCore Trans SLOW	1 × 5 ml
ParaBond Non-Rinse Conditioner + Adhesive A & B	3 × 3 ml
Mixing Tip Short Super Fine	20 pcs

ParaCore SLOW 5ml Intro Kit	6001 3753
ParaCore Dentin SLOW	2 × 5 ml
ParaCore Trans SLOW	2 × 5 ml
ParaBond Non-Rinse Conditioner + Adhesive A & B	3 × 3 ml
Mixing Tip Short Super Fine	20 pcs

ParaCore 5ml White Refill	5886
ParaCore White	2 × 5 ml
Mixing Tip Short Super Fine	20 pcs

ParaCore 5ml Dentin Refill	5887
ParaCore Dentin	2 × 5 ml
Mixing Tip Short Super Fine	20 pcs each

ParaCore 5ml Dentin SLOW Refill	6001 1391
ParaCore Dentin SLOW	2 × 5 ml
Mixing Tip Short Super Fine	20 pcs

ParaCore 5ml Trans SLOW Refill	6001 1392
ParaCore Trans SLOW	2 × 5 ml
Mixing Tip Short Super Fine	20 pcs

Mixing Tip Short Super Fine	6747
Mixing Tip, Ø 1 mm	40 pcs

Mixing Tip Short Fine	6759
Mixing Tip, Ø 1,8 mm	40 pcs

Dual curing core build-up material in 25 ml cartridge.



ParaCore 25ml Dentin Intro Kit	1500
ParaCore Dentin	1 × 25 ml
COLTENE Dispenser 25ml	1 pc
ParaBond Non-Rinse Conditioner + Adhesive A & B	3 × 3 ml
Mixing Tip Yellow / Oral Tip Yellow	20 pcs each

ParaCore 25ml Dentin Kit	1501
ParaCore Dentin	1 × 25 ml
ParaBond Non-Rinse Conditioner + Adhesive A & B	3 × 3 ml
Mixing Tip Yellow / Oral Tip Yellow	20 pcs each

ParaCore 25ml White Intro Kit	1502
ParaCore White	1 × 25 ml
COLTENE Dispenser 25ml	1 pc
ParaBond Non-Rinse Conditioner + Adhesive A & B	3 × 3 ml
Mixing Tip Yellow / Oral Tip Yellow	20 pcs each

ParaCore 25ml White Kit	1503
ParaCore White	1 × 25 ml
ParaBond Non-Rinse Conditioner + Adhesive A & B	3 × 3 ml
Mixing Tip Yellow / Oral Tip Yellow	20 pcs each

ParaCore 25ml Dentin Refill	5853
ParaCore Dentin	1 × 25 ml
Mixing Tip Yellow / Oral Tip Yellow	20 pcs each

ParaCore 25ml White Refill	5854
ParaCore White	1 × 25 ml
Mixing Tip Yellow / Oral Tip Yellow	20 pcs

Mixing Tips Yellow	6550
Mixing Tip Yellow	40 pcs

Oral Tips Yellow	6555
Oral Tip Yellow	100 pcs

COLTENE Dispenser 25 ml	4470
-------------------------	------

ParaBond Adhesive	7486
ParaBond Non-Rinse Conditioner	1 × 3 ml
ParaBond Adhesive A	1 × 3 ml
ParaBond Adhesive B	1 × 3 ml

ParaBond Adhesive A & B Refill	7493
ParaBond Adhesive A	1 × 3 ml
ParaBond Adhesive B	1 × 3 ml

ParaBond Non-Rinse Conditioner Refill	7494
ParaBond Non-Rinse Conditioner	1 × 3 ml

DuoCem

Dual curing resin cement.

Coltène cements have been proving their worth for many years, convincing dentists the world over. DuoCem thus builds on the strengths of these recognised brands, combining them with the benefits of an automix syringe. This combination results in an all-purpose resin cement that offers excellent technical characteristics and can be applied with ease. Thus making DuoCem the ideal cement for routine bonding tasks in the dental practice.



DuoCem Intro Kit	6001 4381
DuoCem Dentin	1 × 5 ml
DuoCem White Opaque	1 × 5 ml
Mixing Tip Short Fine	10 pcs
Mixing Tip Short Super Fine	10 pcs
ParaBond Non-Rinse Conditioner + Adhesive A & B	3 × 3 ml

DuoCem Trial Kit Dentin	6001 4391
DuoCem Dentin	1 × 5 ml
Mixing Tip Short Fine	5 pcs
Mixing Tip Short Super Fine	5 pcs
ParaBond Non-Rinse Conditioner + Adhesive A & B	3 × 3 ml

DuoCem Refill Dentin	6001 4382
DuoCem Dentin	1 × 5 ml
Mixing Tip Short Fine	5 pcs
Mixing Tip Short Super Fine	5 pcs

DuoCem Refill White Opaque	6001 4383
DuoCem White Opaque	1 × 5 ml
Mixing Tip Short Fine	5 pcs
Mixing Tip Short Super Fine	5 pcs



ParaBond Adhesive	7486
ParaBond Non-Rinse Conditioner	1 × 3 ml
ParaBond Adhesive A	1 × 3 ml
ParaBond Adhesive B	1 × 3 ml

ParaBond Adhesive A & B Refill	7493
ParaBond Adhesive A	1 × 3 ml
ParaBond Adhesive B	1 × 3 ml

ParaBond Non-Rinse Conditioner Refill	7494
ParaBond Non-Rinse Conditioner	1 × 3 ml

Mixing Tip Short Super Fine	6747
Mixing Tip, Ø 1 mm	40 pcs

Mixing Tip Short Fine	6759
Mixing Tip, Ø 1,8 mm	40 pcs

Duo Cement Plus

Light- and chemical-cured resin luting cement for cementation of inlays, onlays, adhesive bridges etc.




Duo Cement Plus Kit	7385
Base	2 × 2.5 g
Catalyst	2 × 2.5 g
Accessories	

SELF-ADHESIVE DUAL-CURING RESIN CEMENT



SoloCem[®]

One step reliable cementation

- easy to use and time saving
- high bonding strengths without additional adhesive
- low shrinkage
- antibacterial zinc oxide
- high radiopacity



Intro Kit
REF 6001 3986

SoloCem

Self-adhesive dual curing resin cement.

Secure cementing does not depend on the number of working steps. Self-adhesive SoloCem achieves reliable values – with no additional bonding. This is the result of monomers contained and the composite-like formulation. The convenient auto-mix syringe and the ready-to-use mixing tips offer a simple and time-saving application of SoloCem.



SoloCem Intro Kit	6001 3986
SoloCem Dentin	1 × 5 ml
SoloCem White Opaque	1 × 5 ml
Mixing Tip Short Fine	10 pcs
Mixing Tip Short Super Fine	10 pcs

SoloCem Refill Dentin	6001 4062
SoloCem Dentin	1 × 5 ml
Mixing Tip Short Fine	5 pcs
Mixing Tip Short Super Fine	5 pcs

SoloCem Refill White Opaque	6001 4061
SoloCem White Opaque	1 × 5 ml
Mixing Tip Short Fine	5 pcs
Mixing Tip Short Super Fine	5 pcs

Mixing Tip Short Super Fine	6747
Mixing Tip, Ø 1 mm	40 pcs

Mixing Tip Short Fine	6759
Mixing Tip, Ø 1,8 mm	40 pcs

MAX Pin



Self-threading parapulpal retention pins with shoulder stop and rounded head design for good adaptation. Corrosion-free, titanium alloy. Self-alignment during placement.

Complete Kit 017, blue	M61
Diameter: .017" / .425 mm	
MAX Pin	25 pcs
MAX Drill	2 pcs
Universal hand driver	1 pc

Komplett-Kit 021, purple	M71
Diameter: .021"/.525 mm	
MAX Pin	25 pcs
MAX Drill	2 pcs
Universal hand driver	1 pc

Bulk Kit 021, purple	M73
Diameter: .021"/.525 mm	
MAX Pin	60 pcs



Bulk Kit 017, blue	M63
Diameter: .017"/.425 mm	
MAX Pin	60 pcs

Accessories	
MAX Drills (.017" × .075"/.425 mm × 1.9 mm) blue (depth limiting)	M91
MAX® Drills (.021" × .095"/.525 mm × 2.4 mm) purple (depth limiting)	M92
TMS Regular and Minikin hand wrench (Complete/Two-in-One/Self-Shearing)	B59-3
TMS Minim and Minim 2/1 hand wrench	B59-5
TMS Minuta hand wrench	B59-10
Universal hand driver	L590

TMS LINK / TMS LINK Plus



Self-threading parapulpal retention pins. Alpha-beta titanium alloy or gold plated stainless steel. Shoulder stop and smooth self-alignment during placement. Available as single shear pins and double shear pins. The system is completed with the colour coded KODEX®-drills to ensure an optimal preparation of the canal.



Titanium Alloy

TMS LINK Minikin, red	L811
Single-Shear Pins, Diameter: .017"/.425 mm	
Single-Shear Pins	20 pcs
Kodex drill K91	1 pc
Universal hand driver	1 pc


TMS LINK Plus Minim, silver	EL821-24
Single-Shear Pins Diameter: .021"/.525 mm	
Single-Shear Pins	24 pcs
Kodex drill K92	1 pc
Universal hand driver	1 pc

TMS LINK Plus Minim 2/1, silver	EL841-24
Double-Shear Pins Diameter: .021"/.525 mm	
Double-Shear Pins	24 pcs
Kodex drill K92	1 pc
Universal hand driver	1 pc


TMS LINK Minikin, red	L812
Diameter: .017"/.425 mm (no drills)	
Bulk Kit Single-Shear Pins	50 pcs

TMS LINK Plus Minim, silver	EL822-60
Diameter: .021"/.525 mm (no drills)	
Single-Shear Pins	60 pcs

TMS LINK Plus Minim 2/1 Bulk, silver	EL842-60
Diameter: .021"/.525 mm (no drills)	
Double-Shear Pins	60 pcs



L501



L511

Gold-Plated Stainless Steel

TMS LINK Minuta, pink	L501
Diameter: .0135"/.350 mm	
Single-Shear Pins	20 pcs
Kodex drill K89	1 pc
Universal hand driver	1 pc

TMS LINK Minikin, red	L511
Diameter: .017"/.425 mm	
Single-Shear Pins	20 pcs
Kodex drill K91	1 pc
Universal hand driver	1 pc

TMS LINK Minim, silver	L521
Diameter: .021"/.525 mm	
Single-Shear Pins	20 pcs
Kodex drill K92	1 pc
Universal hand driver	1 pc

TMS LINK Minim 2/1, silver	L531
Diameter: .021"/.525 mm	
Double-Shear Pins	20 pcs
Kodex drill K92	1 pc
Universal hand driver	1 pc

TMS LINK Regular, gold	L541
Diameter: .027"/.675 mm	
Single-Shear Pins	20 pcs
Kodex drill K96	1 pc
Universal hand driver	1 pc

TMS LINK Regular 2/1, gold	L551
Diameter: .027"/.675 mm	
Double-Shear Pins	20 pcs
Kodex drill K96	1 pc
Universal hand driver	1 pc

TMS LINK Minuta Bulk, pink	L502
Diameter: .0135"/.350 mm (no drills)	
Single-Shear Pins	50 pcs


TMS LINK Minikin Bulk, red	L512
Diameter: .017"/.525 mm (no drills)	
Single-Shear Pins	50 pcs

TMS LINK Minim Bulk, silver	L522
Diameter: .021"/.525 mm (no drills)	
Single-Shear Pins	50 pcs


TMS LINK Minim 2/1 Bulk, silver	L532
Diameter: .021"/.525 mm (no drills)	
Double-Shear Pins	50 pcs

TMS LINK Regular Bulk, gold	L542
Diameter: .027"/.675 mm (no drills)	
Single-Shear Pins	50 pcs


TMS LINK Regular 2/1 Bulk, gold	L552
Diameter: .027"/.675 mm (no drills)	
Double-Shear Pins	50 pcs




EL721-24



EL731-24



EL741-24



EL751-24

TMS LINK Plus Minim, silver	EL721-24
Diameter: .021"/.525 mm	
Single-Shear Pins	24 pcs
Kodex drill K92	1 pc
Universal hand driver	1 pc

TMS LINK Plus Minim 2/1, silver	EL731-24
Diameter: .021"/.525 mm	
Double-Shear Pins	24 pcs
Kodex drill K92	2 pcs
Universal hand driver	1 pc

TMS LINK Plus Regular, gold	EL741-24
Diameter: .027"/.675 mm	
Single-Shear Pins	24 pcs
Kodex drill K96	1 pc
Universal hand driver	1 pc

TMS LINK Plus Regular 2/1, gold	EL751-24
Diameter: .027"/.675 mm	
Double-Shear Pins	24 pcs
Kodex drill K96	2 pcs
Universal hand driver	1 pc

TMS LINK Plus Minim Bulk, silver	EL722-60
Diameter: .021"/.525 mm (no drills)	
Single-Shear Pins	60 pcs

TMS LINK Plus Minim 2/1 Bulk, silver	L732-60
Diameter: .021"/.525 mm (no drills)	
Double-Shear Pins	60 pcs

TMS LINK Plus Regular Bulk, gold	EL742-60
Diameter: .027"/.675 mm (no drills)	
Single-Shear Pins	60 pcs

TMS LINK Plus Regular 2/1 Bulk, gold	EL752-60
Diameter: .027"/.675 mm (no drills)	
Double-Shear Pins	60 pcs

KODEX-Drills



Minuta pink	K89
.0135" x .051"/.350 x 1.3 mm, depth-limiting	
	6 pcs

Minikin red	K90
.017" x .059"/.425 x 1.5 mm, short shank depth-limiting	
	6 pcs

Minim silver	K92
.021" x .079"/.525 x 2 mm, depth-limiting	
	6 pcs

Regular gold	K96
.027" x .079"/.675 x 2 mm, depth-limiting	
	6 pcs

Minikin red	K91
.017" x .059"/.425 x 1.5 mm, long shank, depth-limiting	
	6 pcs

Minim silver	K93
.021" x .157"/.525 x 4 mm, without depth-limitation	
	6 pcs

Regular gold	K97
.027" x .197"/.675 x 5 mm, without depth-limitation	
	6 pcs

Universal hand driver	L590
	1 pc

TMS



TMS Regular Complete Kit	T01
Diameter (.027"/.675 mm)	
Self-Threading Pins	20 pcs
Kodex Drill K97	1 pc
Hand Wrench	1 pc

TMS Regular Refill Kit	T02
Diameter (.027"/.675 mm)	
Self-Threading Pins	50 pcs



TMS Regular Bulk Kit	T03
Diameter (.027"/.675 mm)	
Self-Threading Pins	100 pcs

TMS Minim Refill Kit	T05
Diameter (.021"/.525 mm)	
Self-Threading Pins	50 pcs



TMS Minim Complete Kit	T04
Diameter (.021"/.525 mm)	
Self-Threading Pins	20 pcs
Kodex Drill K93	1 pc
Hand Wrench	1 pc

TMS Minim Bulk Kit	T06
Diameter (.021"/.525 mm)	
Self-Threading Pins	100 pcs

TMS Regular Two-in-One Kit	T07
Diameter (.027"/.675 mm)	
Self-Threading Pins	2 x 20 pcs
Kodex Drill K96	1 pc
Hand Wrench	1 pc

TMS Regular Two-in-One Refill Kit	T08
Diameter (.027"/.675 mm)	
Self-Threading Pins	2 x 50 pcs



TMS Regular Two-in-One Bulk Kit	T09
Diameter (.027"/.675 mm)	
Self-Threading Pins	2 x 100 pcs

TMS Minim Two-in-One Refill Kit	T11
Diameter (.021"/.525 mm)	
Self-Threading Pins	2 x 50 pcs

TMS Minim Two-in-One Kit	T10
Diameter (.021"/.525 mm)	
Self-Threading Pins	2 x 20 pcs
Kodex Drill K92	1 pc
Hand Wrench	1 pc

TMS Minim Two-in-One Bulk Kit	T12
Diameter (.021"/.525 mm)	
Self-Threading Pins	2 x 100 pcs



TMS Regular Self-Shearing Complete Kit		T13
Diameter (.027"/.675 mm)		
Self-Threading Pins	20 pcs	
Kodex Drill K97	1 pc	
Hand Wrench	1 pc	

TMS Regular Self-Shearing Refill Kit		T14
Diameter (.027"/.675 mm)		
Self-Threading Pins	50 pcs	



TMS Regular Self-Shearing Bulk Kit		T15
Diameter (.027"/.675 mm)		
Self-Threading Pins	100 pcs	



TMS Minikin Complete Kit		T107
Diameter (.017"/.425 mm)		
Self-Threading Pins	20 pcs	
Kodex Drill K91	1 pc	
Hand Wrench	1 pc	





TMS Minikin Refill Kit		T108
Diameter (.017"/.425 mm)		
Self-Threading Pins	50 pcs	



TMS Minikin Bulk Kit		T109
Diameter (.017"/.425 mm)		
Self-Threading Pins	100 pcs	

TMS Minuta Complete Kit		T111
Diameter (.0135"/.350 mm)		
Self-Threading Pins	20 pcs	
Kodex Drill K89	1 pc	
Hand Wrench	1 pc	

TMS Minuta Refill Kit		T112
Diameter (.0135"/.350 mm)		
Self-Threading Pins	50 pcs	



CH **Coltène/Whaledent AG**
 Feldwiesenstr. 20
 9450 Altstätten / Switzerland
 T + 41 71 757 53 00
 F + 41 71 757 53 01
 info.ch@coltene.com

CN **Coltene/Whaledent China Trading Co. Ltd.**
 Rm. 1103/1104, Bldg. B2, Wu Dong Plaza
 #9 Chegongzhuang Street, West District
 Beijing 100044 / P.R. China
 T + 86 10 6212 3251
 F + 86 10 6212 3252
 info.cn@coltene.com

IT **Coltène Italia S.r.l.**
 c/o Polo Tecnologico della Brianza
 Via Alessandro Volta, 94 - Cond. Bianchi
 20832 Desio (MB) / Italy
 T + 39 0362 62 67 32
 F + 39 0362 63 25 80
 info.it@coltene.com

DE **Coltène/Whaledent GmbH + Co. KG**
 Raiffeisenstraße 30
 89129 Langenau / Germany
 T + 49 7345 805 0
 F + 49 7345 805 201
 info.de@coltene.com

FR **Coltène/Whaledent s.a.r.l.**
 75 Boulevard Marie et Alexandre Oyon
 72058 Le Mans Cedex 2 / France
 T + 33 2 43 39 30 30
 F + 33 2 43 39 30 40
 info.fr@coltene.com

UK **Coltène/Whaledent Ltd.**
 The President Suite
 Kendal House, Burgess Hill
 Victoria Way
 West Sussex, RH15 9NF / UK
 T + 44 1444 235486
 F + 44 1444 870640
 info.uk@coltene.com

US **Coltène/Whaledent Inc.**
 235 Ascot Parkway
 Cuyahoga Falls, Ohio 44223 / USA
 T + 1 330 916 8800
 F + 1 330 916 7077
 info.us@coltene.com

IN **Coltène/Whaledent Private Ltd.**
 106/107, Hallmark, Vasant Oscar
 L.B.S Marg, Mulund West
 Mumbai 400080, Maharashtra / India
 T + 91 22 25923626
 F + 91 22 25623165
 info.in@coltene.com

BR **Vigodent SA Indústria e Comércio**
 Rua Pesqueira 26 Bonsucesso
 Rio de Janeiro CEP 21041-150 / Brazil
 T + 5521 3865 5600
 F + 5521 2280 8411
 sac@vigodent.com.br

© 2015 Coltène/Whaledent AG – www.coltene.com



COLTENE

