



Whisper™

Pure Sapphire

التقويم الكريستالي الشفاف

We're proud to offer Intrigue Whisper™, the highest quality aesthetic bracket on the market. Polycrystalline material makes Intrigue Whisper the only crystal-clear braces available. They're not only discreet, they're virtually invisible.

- ◆ Polycrystalline 99.99% Pure Transparency
- ◆ Invisible Clear
- ◆ Improved staining resistance
- ◆ Torque in base
- ◆ Maximum tie wing strength
- ◆ High strength by maintaining flexibility for improved shock resistance
- ◆ Patented mechanical bonding base for strong bonds yet easy de-bonding
- ◆ Available Roth & MBT Technique .018/.022 Slot
- ◆ High tech material increased strength of 10%
- ◆ Low-friction Polyurethane for superior strength
- ◆ Abrasion & friction resistance
- ◆ Easy to ligate
- ◆ Smooth & Comfortable



*Your choice for sensible
orthodontic solution*

